



Plasmid name : pQCXIPpuro H3.3-K4M

pQCXIP Vector Information

PT3669-5
Cat. No. 631516



Fragment Insert (indicate source, restriction ends, whether sequenced, if PCR derived) : _____

PCR : H3.3-XhoI-K4M(5') + H3.3-EcoR1(3'); pQCXIPpuro-H3.3-WT-Flag-HA as a template

Antibiotic resistance : Amp

Note : Double digestion; XhoI/EcoR1 (7.2kb+0.5kb)
Confirmed by sequencing and Western blot

References : _____

Constructed by Young Hoon Jang Date : 8-25-2015