

pUC57-PB-PGK-tetR-3xFLAG-HA-2A-Puro

5518 bp

PB left 5' IR 446..480

PB left ID 481..763

hPGK 769..1287

PGK 764..1322

tetR-NLS_optimal 1323..2066

3xFlag 2073..2141

HA 2148..2174

GSG for T2A 2196..2204

T2A 2199..2261

T2A core 2205..2258

Puro-optimal 2262..2861

PB right 3' IR 3159..3221

PB right ID 2982..3158

SV40 polyA 2862..2981

