


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLJM1-Tmem192-mRFP-3xHA with the following enzymes: EcoRI-HF, EcoRV-HF, and NheI-HF.

Plasmid 134631:
pLJM1-Tmem192-mRFP-3xHA



DNA digested with NheI-HF/EcoRI-HF and EcoRV-HF/EcoRI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of SYBR Safe.

Ran for 1.5 hours at 100 volts.