


Please submit an electronic copy of this document to the DNA librarian. Provide an aliquot of $\geq 1\mu\text{g}$ (low copy) or $\geq 5\mu\text{g}$ (high copy) DNA at $100\text{ng}/\mu\text{L}$ in a 1.5mL tube labeled with the construct name and date. (Note: low stocks should be resupplied by the person using the DNA).

DNA Library #	Name of the Construct (self explanatory, no code please!)	Box number/ position	Date Deposited	Source (Person who created it or the person/lab that sent it to us)	Person who deposited in the library (a name from our lab)
---------------	---	----------------------	----------------	---	---



Proteins or fragments included:

Purpose of the Plasmid:

Concentration:

Backbone Vector:

Antibiotic Selection Marker:

Any obligations? (any requirements for disclosure, etc)

Species:

Add'l Notes or Comments: (include mutations and functional consequences)

Viral construct?

What type?

Sequenced?

Sequence file sent (.gb and .txt)?

Cloning Site Sequence:

