


CdA (Addgene) T7 RNA Polymerase - Hyperactive AID - Fragment 1 - CdA (Addgene) Copy - Fragment 1 (concatenated sequence 1) (concatenated sequence 3)  
 Linker+UGI (concatenated sequence 2)



pcDNA3.1 Hyperactive AID - T7 RNA Polymerase - UGI Copy - Fragment 1 - pcDNA3.1 IRES+EGFP - Fragment 1 (concatenated sequence 1) (concatenated sequence 3)  
 CdA (Addgene) T7 RNA Polymerase Copy - Fragment 1 - CdA (Addgene) Copy - Fragment 1 (concatenated sequence 1) (concatenated sequence 1)