


This gel depicts the NEB 1kb ladder and a digestion of the plasmid 14715 with the following enzymes: XmaI single digest, AgeI-HF single digest, AgeI-HF + XmaI double digest, EcoRI-HF + BamHI-HF double digest

Plasmid 14715
Lentiviral-TOP-dGFP-reporter



DNA digested with XmaI single digest, AgeI-HF single digest, AgeI-HF + XmaI double digest, EcoRI-HF + BamHI-HF double digest

for **1.5 hours**

Gel: 0.8% TAE stained with 100 ml 1X TAE and 5uL of ethidium bromide

Ran for 2 hours at **47** volts