


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLVX-EF1alpha-SARS-CoV-2-orf3a-2xStrep-IRES-Puro with the following enzymes: EcoRI-HF and MluI-HF.

Plasmid 141383:
pLVX-EF1alpha-SARS-CoV-2-orf3a-
2xStrep-IRES-Puro



DNA digested with EcoRI-HF and MluI-HF single digests and EcoRI-HF/ MluI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.