

Glycerol stock # 987

DNA: 4°C Kalirin 2 Box #59  
-70°C Ping # Box #52

INVESTIGATOR: Yanping Wang

DATE: 12/13/2012


pEAK-rat cKal7

**E-pDNA INVENTORY SHEET**

Clone Designation: pEAK-rKal7c Cloned DNA Source: PCR  
Rat brain cortex cDNA as a template  
Description: VECTOR: pEAK-Hismyc-rKal7a(93nt) Insertion Site: NcoI/NruI  
INSERT SIZE (Kb): 145 nt  
Bacterial Host: Top10 Strain Designation: \_\_\_\_\_  
Growth Medium: LB  
Antibiotic Selection: Amp

**Restriction Map**

Genbank AF230644.1



[J Mol Neurosci](#). 2004;22(1-2):51-62.

**Kalirin expression is regulated by multiple promoters.**

[McPherson CE](#)<sup>1</sup>, [Eipper BA](#), [Mains RE](#).