


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid pSLQ5429_pUC_hU6-crScaffold_EF1a-BFP with the following enzymes: EcoRI-HF, Xho I

Plasmid 155306:
pSLQ5429_pUC_hU6-crScaffold_EF1a-BFP



DNA digested with EcoRI-HF single digest and EcoRI-HF/Xho I double digest for **2 hours**.

Gel: 0.8% agarose stained with SYBR Safe.

Ran for 1.25 hours at 100 volts.