


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pTol2 hsp70l loxP-DsRed-loxP-nCas9n-GFP with the following enzymes: SacI-HF and EcoRV-HF.

Plasmid 158962:

pTol2 hsp70l loxP-DsRed-loxP-nCas9n-GFP



DNA digested with a SacI-HF/EcoRV-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.