


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pCFJ2334 - [1-2] ENTR - Fluor - ce-GFP(PATCs(900bp), NLS, no_atg, no stop) with the following enzyme: BamHI.

Plasmid 159850:

pCFJ2334 - [1-2] ENTR - Fluor - ce-GFP(PATCs(900bp), NLS, no_atg, no stop)



DNA digested with a BamHI single digest for **2 hours**.
Gel: 0.8% agarose stained with 5 uL of ethidium bromide.
Ran for 1.5 hours at 100 volts.