



pUltra-Plus-H2B-EGFP



Nuclear expression of the H2B-EGFP gene in BHK-21 cells at 24 hours post-transfection with the plasmid pUltra-Plus-H2B-EGFP. Total magnification of 200X (A) and 600X (B). Fluorescence live-cell images were taken with a Lionheart FX Automated Microscope (BioTek) using the GFP Filter / 465 nm LED cube set. Images by Jorge L. Arias-Arias (jorgeluis.arias@ucr.ac.cr), Centro de Investigación en Enfermedades Tropicales (CIET), Universidad de Costa Rica (<http://www.ciet.ucr.ac.cr/>).