


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pFlox\_N-p300\_Neo\_FKBPF36V\_mScarlet with the following enzymes:HindIII-HF and Scal-HF

**Plasmid 163886:**  
pFlox\_N-p300\_Neo\_FKBPF36V\_mScarlet



DNA digested with HindIII-HF/ Scal-HF double digests for **2 hours**.  
Gel: 0.8% agarose stained with 5 uL of ethidium bromide.  
Ran for 1.5 hours at 100 volts.