


pMSCVhyg_T7_CA-Chrebp plasmid



Insert: PCR primers: XhoI_mCA-Chrebp_FW + HpaI_mCA-Chrebp_RV;
N-terminal 3xT7 as a template

Vector: pMSCVhyg_N_3xT7+2xStre → XhoI + HpaI

Antibiotic resistance: Amp

Note: XhoI + HpaI + 3xT7 + 2xStre + Chrebp → 1.8 kb

Confirmed by sequencing + Western blot

References: pMSCVhyg vector from Addgene
Constructed by Nhien Tran, 01-29-2019