


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid pAAV2/1-CAG-Flex-FlpO-WPRE with the following enzymes: HindIII-HF and PvuI-HF.

Plasmid 171690
pAAV2/1-CAG-Flex-FlpO-WPRE



DNA digested with HindIII-HF and PvuI-HF single digests and HindIII-HF/PvuI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.