


Gene name- pCS2+ Wbp2nl

Vector- PCS2+

Purpose- mRNA expression



Vector size: 4.1 Kb

Insert size- about 948bp

Selection antibiotic: AMP

Cloning Sites: Cla1/Xho1

Tag-None

Constructed by: Diana Johnson

Date:

GenBank Accession Number: BC082812.1

Reference:

ANTISENSE: Restriction enzyme: Cla1

Polymerase: T3

SENSE: Restriction enzyme: Not1

Polymerase: SP6

COMMENTS: -66 5' UTR, full ORF and no 3' UTR. Morpholino sensitive.