


Gene – pCS2+ 5'MT-Wbp2nl  
Vector- PCS2+MT  
Purpose- mRNA expression



Vector size: 4.3 Kb  
Insert size- about 882bp  
Selection antibiotic: AMP  
Cloning Sites: EcoR1/Xho1  
Tag- Myc

Constructed by: Himani Datta Majumdar  
Date:  
GenBank Accession Number:

Reference:

ANTISENSE:	Restriction enzyme: EcoR1	Polymerase: T3
SENSE:	Restriction enzyme: Not1	Polymerase: SP6

COMMENTS: myc tag is at N-terminus, insert does not contain endogenous start site.  
Morpholino insensitive. No 5' & 3' UTRs.