


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid AAV-FLEX-rev-ChR2-tdtomato with the following enzymes: EagI-HF, NotI-HF, Scal-HF

**Plasmid 18917:**  
AAV-FLEX-rev-ChR2-tdtomato



DNA digested with EagI-HF single digest and NotI-HF/Scal-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5  $\mu$ L of ethidium bromide.

Ran for 1.5 hours at 100 volts.