


# Plasmid name: pCW57.1-Cdx2-3xT7

Created with SnapGene®



Template	pLV-TetO-Cdx2 (Addgene #70697) for Cdx2 CDS; pMSCVhyg-3xT7 (Plasmid #171028) for 3xT7
Insert approach	PCR by MluI-Cdx2 (5') + BamHI-Cdx2 (3') and BamHI-3xT7 (5') + HpaI-3xT7 (3')
Vector	pCW57.1_Gata6 (Addgene #73537), replaced Gata6 CDS with multiple cloning sites by mutagenesis
Antibiotic resistance	Amp
Note	Confirmed by sequencing
References	