


Ligation of 2 sequences.



AAAA-0245-pLV-RnSyt1-sgRNA1-dCas9-KRAB-TagBFP
15,009 bp

DNA Type: Synthetic DNA

Unspecified
Dam⁻ Dcm⁻ EcoKI⁻

Description: Ligation of 2 sequences.

Created: 12 Feb 2021

Last Modified: 12 Feb 2021

Accession Number: urn.local...bc-crn8qv7

Code Number:

Sequence Author:

Comments:

References:

Embedded Files: