


This gel depicts the NEB 1kb extend ladder and a digestion of the plasmid piggyBac-rtTA (4th_Gen)-sfGFP-IRES-NGN2-puro (VK_1122) with the following enzyme(s): KpnI-HF and PvuI-HF.

Plasmid 209079:

piggyBac-rtTA (4th_Gen)-sfGFP-IRES-NGN2-puro (VK_1122)



DNA digested with KpnI-HF and PvuI-HF single digests for **2 hours**.
Gel: 0.8% agarose stained with 5 uL of SYBR Safe.
Ran for 1.5 hours at 100 volts.