


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pBabe-TetCMV-puro-mCherry-PA-Rac1 with the following enzymes: EcoRI-HF and XhoI.

**Plasmid 22035:**  
pBabe-TetCMV-puro-mCherry-PA-Rac1



DNA digested with EcoRI-HF and XhoI single digests and EcoRI-HF/XhoI double digest for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.