





K4-PB-TO_APEX-V5-NES (14104 bp)




















taataaattcaacaaacaattttatttattgtttattttattttaaaaaaacaacaaactcaaaatttcttctataaagtaacaaaacttttatgagggacagcccc
attatttaagtgtgtttgtaataataacaaataataataattttttttgtttttgagttttaagaagatatttcattgttttgaaaatactccctgtcggggg
gttttgaaaatactccctgtcgg











tatcagtgatagagaacgtatctacagtttactccctatcagtgatagagaacgtatatccagtttactccctatcagtgatagagaacgtataagctttaggcgtg
atagtcactatctcttgcatagatgtcaaatgagggatagtcactatctcttgcataaggtcaaatgagggatagtcactatctcttgcataatcgaaatccgcac



3,320 3,330 3,340 3,350 3,360 3,370 3,380 3,390 3,400 3,410 3,420






tacggtgggcgccctataaaagcagagctcggttagtgaaccgtcagatcgccctggagcaattccacaacacttttgtcttataccaactttccgtaccacttcttac
atgccaccgcggatattttcgctcgcagcaaatcacttggcagtcctagcggacctcgtaaggtgtgtgaaaacagaatatggttgaaaggcatggtgaaggatg



3,430 3,440 3,450 3,460 3,470 3,480 3,490 3,500 3,510 3,520 3,530



cctcgtaaaggatccgccgccaccatgGGAAGTCTTACC
cctcgtaaaggatccgccgccaccatgGGAAGTCTTACCAACTGTGAGTGTCTGATTACCAGGACGCCGTTGAGAAGGCGAAGAAGAAGCTCAGAGGCTTCATCGC
ggagcatttcctaggcggcggtggtacCCTTTCAGAATGGTTGACACTACGACTAATGGTCTCGGCAACTCTTCCGTTCTTCTTCGAGTCTCCGAAGTAGCG



3,540 3,550 3,560 3,570 3,580 3,590 3,600 3,610 3,620 3,630

Abasi Abasi Abasi Abasi Abasi

TGAGAAGAGATGCGCTCCTCTAATGCTCCGTTTGGCATTCCACTCTGCTGGAACCTTTGACAAGGGCACGAAGACCGGTGGACCCTTCGGAACCATCAAGCACCTG
ACTCTTCTCTACGCGAGGAGATTACGAGGCAAACCGTAAGGTGAGACGACCTTGAAACTGTTCCCGTGCTTCTGGCCACCTGGGAAGCCTTGGTAGTTCTGTTGGGAC

GGAGATTACGAGGCAAACCG

28 30 32 34 36 38 40 42 44 46 48 50 52 54 56 58 60 62

E K R C A P L M L R L A F H S A G T F D K G T K T G G P F G T I K H P

3559-4386

source

synthesized

SeqVer (WT AAV R)

SeqVer (TRE 3G F)

source

SeqVer (SV40 F)

pcDNA3 APEX-NES

APEX2

3,640 3,650 3,660 3,670 3,680 3,690 3,700 3,710 3,720 3,730 3,740

Abasi AfeI Abasi Abasi Abasi Abasi AciI Abasi Abasi

CCGAAGTGGCTCACAGCGCTAACAAACGGTCTTGACATCGCTGTTAGGCTTTTGGAGCCACTCAAGGCGGAGTTCCCTATTTTGGAGCTACGCCGATTCTACCAGTTG
GGCTTGACCGAGTGTGCGGATTGTTGCCAGAACTGTAGCGACAATCCGAAAACCTCGGTGAGTTCGCCCTCAAGGGATAAACTCGATGCGGCTAAAGATGGTCAAC

64 66 68 70 72 74 76 78 80 82 84 86 88 90 92 94 96 98

A E L A H S A N N G L D I A V R L L E P L K A E F P I L S Y A D F Y Q L

3559-4386

source

synthesized

SeqVer (WT AAV R)

SeqVer (TRE 3G F)

source

SeqVer (SV40 F)

pcDNA3 APEX-NES

APEX2

3,750 3,760 3,770 3,780 3,790 3,800 3,810 3,820 3,830 3,840 3,850

