



Plasmid name: pCX-Flag-P/CAF

Plasmid size: 7.00 kb

Constructed by: Xiang-Jiao Yang and Yoshihiro Nakatani

Construction date: April 1995; restriction sites are rough estimates.

Comments/References: The P/CAF cassette flanked by EcoRI just upstream from the Flag sequence and BamHI downstream from the polyadenylation signal was inserted between EcoRI and BglIII sites of pCX. Yang et al. (1996) Nature 328, 319; Minucci et al. (1994) MCB 14, 360.