


This gel depicts the NEB 1kb ladder and a digestion of the plasmid 26681 with the following enzymes: SmaI single digest, EcoRI-HF Single digest, and SmaI and EcoRI-HF double digest

**Plasmid 26681**  
pCMVTNT-T7-KPNA5



DNA digested with SmaI single digest, EcoRI-HF Single digest, SmaI and EcoRI-HF double digest

**2 hours**

Gel: 0.8% TAE with 3  $\mu$ L of ethidium bromide in the gel.

Ran for 4 hours at **47** volts