


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pDD162 (Pefft-3::Cas9 + Empty sgRNA) with the following enzymes: SpeI-HF and NcoI-HF.

Plasmid 47549:
pDD162 (Pefft-3::Cas9 + Empty sgRNA)



DNA digested with SpeI-HF and NcoI-HF double digest for **2 hours**.
Gel: 0.8% agarose stained with 6 uL of ethidium bromide.
Ran for 1.5 hours at 100 volts.