

This gel depicts the NEB 1kb ladder and a digestion of the plasmid 48742 with the following enzymes: Sall-HF single digest, XhoI single digest, BamHI-HF single digest, and Sall-HF and XhoI double digest

Plasmid 48742

pGL2

Sall-HF and
XhoI double
48742


BamHI-HF 48742

XhoI 48742

Sall-HF 48742

Uncut 48742

1kb Ladder



DNA digested with Sall-HF single digest, XhoI single digest, BamHI-HF single digest, Sall-HF and XhoI double digest

2 hours

Gel: 0.8% TAE with 3 uL of ethidium bromide in the gel.

Ran for 3 hours at **47** volts