


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLKO.1-puro U6 sgRNA BfuAI large stuffer with the following enzymes: KpnI-HF and KpnI-HF/EcoRI-HF.

Plasmid 52628:
pLKO.1-puro U6 sgRNA BfuAI large stuffer



DNA digested with KpnI-HF single digest and KpnI-HF/EcoRI-HF double digest for **2 hours.**

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 2 hours at 80 volts.