


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pZac2.1 gfaABC1D-lck-GCaMP6f with the following enzymes: PvuI, KpnI HF and NotI HF

**Plasmid 52924:**  
pZac2.1 gfaABC1D-lck-GCaMP6f



DNA digested with PvuI, KpnI HF and NotI HF single digests for **2 hours**  
Gel: .8% TAE stained with 150 ml 1X TAE and 3uL of ethidium bromide  
Ran for 1.5 hours at 100 volts