

Reporter H4

Box I reporter G6

Plasmid Name: pIhhP-1.3k-luc Size: 7.9 k

Made by: Weiguang Wang Date: 7-17-06

Vector: P₄-luc Size: 6.6 k

Insert: IhhP 1.3 k Size: 1.3 k

Cloning sites: HpaI, HindIII Selection: Amp^r

Notes: 1.3 k Ihh promoter fragment is cut from pIhhP2.8 k by HpaI and HindIII, then insert into P₄luc in SmaI and HindIII cutting site.

