

Reporter Hg

Plasmid Name: pLUC-Ihhp742 Size: 7.39 k
 Made by: Weiguang Wang Date: 2/24/2006
 Vector: p4 LUC Size: 6.6 k
 Insert: Ihh promoter 742 bp Size: 742 bp
 Cloning sites: XhoI & NotI Selection: Amp⁺
 Notes: Ihhp742 bp promoter fragment (PCR Primer, P₂ = CTG-
 AGA AAG GGA ATG TTG CC, P₄ = GCG TGC TGT CCCC CT CGG CG)
 was inserted into p4LUC in XhoI and NotI sites.

