


This gel depicts the NEB 1kb ladder and a digestion of the plasmid mEmerald-CLIP170-N-18 with the following enzymes: BamHI-HF, NheI-HF, and Bgl II.

Plasmid 54044:
mEmerald-CLIP170-N-18



DNA digested with BamHI-HF, NheI-HF, and Bgl II single digests and BamHI-HF/NheI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.