


NAME OF CONSTRUCT: pcDNA3-WT-VLCAD

Plasmid backbone: pcDNA3

Antibiotic Resistance: Amp and Neo

NM\_000018.3 Homo sapiens VLCAD



The insert is amplified by using following PCR primers

VLCAD (EcoRI)\_up\_primer: GAATTCATGCAGGCGGCTCGGATGGC

VLCAD (NotI)\_down\_primer: GCGGCCGCTCTGAAGCCAAGTGGGTTGCTG

(template: cDNA, a human fetus brain cDNA library)