

Category	Mammalian expression	File No.	
Plasmid	pRK-FLAG-TRAF4	A.K.A.	FLAG-TRAF4 (FL)
Description	FLAG tagged mouse TRAF4		
Reference	Wang et al, JBC 2013, 288, 21784-92		
Insert	TRAF4 (1-470)	Tag (end)	FLAG(N)
Vector	pRK2F	Vector sites	BamHI/Sall
Expression	Mammalian	Selection	Amp
Constructed by	Chaoyang Jin	Date	9/20/2006
Method Construction	PCR mouse TRAF4 NM_009423 (mage Clone 3598556, GenBank BC020455) with Primers: BamHI-mTraf4-F, and Sall-mTraf4-r; digest PCR produce with BamHI/Sall, ligate into pRK2F		

Sequence information FLAG-EcoRI-BamHI-TRAF4 (FL)-stop-Sall-HindIII

(BamHI) ggatcc ATG CCC GGCAGC TGA gtcgac (Sall)
TRAF4 (FL)

