

Category	Mammalian expression	File No.	
Plasmid	pEGFP-TRAF4(K190R)	A.K.A.	EGFP-TRAF4 (KR)
Description	GFP tagged mouse TRAF4 K190R		
Reference	Wang et al, JBC 2013, 288, 21784-92		
Insert	TRAF4 (K190R)	Tag (end)	GFP(N)
Vector	pEGFP-C1	Vector sites	BglIII/Sall
Expression	Mammalian	Selection	Kan
Constructed by	Ying Zhang	Date	3/29/2012
Method Construction	BamHI/Sall Fragment from pRK-FLAG-TRAF4 (K190R) into BglIII/Sall Fragment of pEGFP-C1. as a result the BglIII site is destroyed.		

Sequence information 20120329: YiZ_EGFP-Traf4KR

GFP-TRAF4(K190R)-stop-SalI

agatcc ATG CCC GGC ...CGG...AGC TGA gtcgac (SalI)

TRAF4 (K190R)

