


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLL7.0: hITSN1(1159-1509)-tgRFPt-SSPB R73Q with the following enzymes: NotI-HF.

**Plasmid 60420:**  
pLL7.0: hITSN1(1159-1509)-tgRFPt-  
SSPB R73Q



DNA digested with NotI-HF single digest for **2 hours**.  
Gel: 0.8% agarose stained with 6 uL of ethidium bromide.  
Ran for 1.5 hours at 100 volts.