


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pX601-AAV-CMV::NLS-SaCas9-NLS-3xHA-bGHpA;U6::Bsal-sgRNA with the following enzymes: SacI and StuI.

**Plasmid 61591:**

pX601-AAV-CMV::NLS-SaCas9-NLS-  
3xHA-bGHpA;U6::Bsal-sgRNA



DNA digested with SacI and StuI single digests for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 2 hours at 80 volts.