


* fully sequenced *

<p>RHD Abbreviated Protein Name</p>	<p>570 Tube Number</p>	<p>pTR1EX - Ritoc biosensor WT Full Name</p>	<p>● DNA ○ RNA ○ Cells (strain)</p>	<p>6.5.08 prep Date created or obtained</p>	<p>H.B8 Location</p>
-----------------------------------------	----------------------------	--------------------------------------------------	----------------------------------------------	-----------------------------------------------------	--------------------------



Purpose mammalian expression

Parental Plasmid pTR1EX (Novogen)

Antibiotic Selection Marker bact = CARB / Amp
mamm = None

Concentration _____

Source Louis HODGSON / JON ZAWISTOWSKI

Comments for transient transfection of
Ritoc biosensor

* T153M change is present in
mVenus for all Ritoc biosensor
derivatives *

Cloning Site Sequence :

