


* fully sequenced *

<p>RHO Abbreviated Protein Name</p>	<p>571 Tube Number</p>	<p>pTRIEX - Ritoc biosensor WT Q63L Full Name</p>	<p>● DNA ○ RNA ○ Cells (strain)</p>	<p>6.5.08 prep Date created or obtained</p>	<p>H.B9 Location</p>
---	----------------------------	--	--	---	--------------------------



Purpose mammalian expression

Parental Plasmid pTRIEX (Novogen)

Antibiotic Selection Marker bact = CARB / Amp
mamm = None

Concentration _____

Source Louis HOGSON / JON ZAWISTOWSKI

Comments for transient transfection of
Ritoc biosensor

* T153M change is present in
mVenus for all Ritoc biosensor
derivatives *

Cloning Site Sequence :

