


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLV-CMV-LoxP-DsRed-LoxP-eGFP with the following enzymes: EcoRV-HF and BamHI-HF.

Plasmid 65726:
pLV-CMV-LoxP-DsRed-LoxP-eGFP



DNA digested with EcoRV-HF and BamHI-HF single digests for **2 hours**.
Gel: 0.8% agarose stained with 6 uL of ethidium bromide.
Ran for 1.5 hours at 100 volts.