


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pKLV2-U6gRNA5(gGFP)-PGKBFP2AGFP-W with the following enzymes: KpnI-HF and XhoI.

Plasmid 67980:
pKLV2-U6gRNA5(gGFP)-
PGKBFP2AGFP-W



DNA digested with KpnI-HF and XhoI single digests and KpnI-HF/ XhoI double digest for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.