


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid 5xgRNA: PTEN exon 5, p53 exon 8, SMAD4 exon 2, p53 exon 7, SMAD4 exon 9 probasin\_mTQ2\_FlpO with the following enzymes: AgeI-HF, SpeI-HF

### Plasmid 68357:

5xgRNA: PTEN exon 5, p53 exon 8, SMAD4 exon 2, p53 exon 7, SMAD4 exon 9 probasin\_mTQ2\_FlpO



DNA digested with AgeI-HF and SpeI-HF double digest for **2 hours**.  
Gel: 0.8% agarose stained with 5 uL of ethidium bromide.  
Ran for 1.5 hours at 100 volts.