


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pAAV-EF1a-C-CreintG with the following enzymes: BamHI-HF and EcoRI-HF.

**Plasmid 69571:  
pAAV-EF1a-C-CreintG**



DNA digested with BamHI-HF and EcoRI-HF single digests and BamHI-HF/EcoRI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5  $\mu$ L of ethidium bromide.

Ran for 1.5 hours at 100 volts.