


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLV hU6-sgRNA hUbC-dCas9-KRAB-T2a-Puro with the following enzymes: NotI-HF and KpnI-HF.

**Plasmid 71236:**  
pLV hU6-sgRNA hUbC-dCas9-KRAB-  
T2a-Puro



DNA digested with NotI-HF/ KpnI-HF double digest for **2 hours**.  
Gel: 0.8% agarose stained with 6 uL of ethidium bromide.  
Ran for 1.5 hours at 100 volts.