

Plasmid name: hHSF1-myc his PRD Size: 7,1 kb

Vector: pcDNA 3.1(-)A/myc-his Size: ~5,5 kb

Insert: human HSF1 Size: ~1,5 kb

Cloning site: EcoRI / HindIII

Epitope tag: myc his Resistance: Ampicillin/Neomycin

Expression: Mammalian: X Bacterial: _____

Glycerol stock: _____

Note: _____

hHSF1-Myc-His

Vector: pcDNA3.1/Myc-His(-)A.

Size: ~7,1 kb

Construction: Vector digested with EcoRI and HindIII and HSF1 ligated to the vector.

1587 nt of hHSF1 coding region

ATG present, stop codon mutated

