


This gel depicts the NEB 1kb ladder and a digestion of the plasmid AAV Synapsin Intron H2B Gcamp 6s wpre Pzac2.1 with the following enzymes: Scal-HF and HindIII-HF.

Plasmid 74150:

AAV Synapsin Intron H2B Gcamp 6s wpre Pzac2.1



DNA digested with a Scal-HF/HindIII-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.