

PLASMID INFORMATION SHEET

PLASMID NUMBER: pp395

PLASMID NAME: FLAG-hPAX3

USED FOR: overexpression (cell culture)

HOST: XL-Blue

SELECTABLE MARKERS: Kan^R

PLASMID HISTORY: 1) Source: MSA; hPAX3 EST clone (Open Biosystems®)

2) Construction: PCR w/cut sites added

a) Insert: pp75 (BamHI/EcoRI) (hPAX3 EST clone) 484 aa

b) Vector: pp191 (Tag2B) (BamHI/EcoRI)

3) Comments 484aa ORF (longest hPAX3 isoform)

