


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pAAV-CAG-FLEX-rc [Chronos-tdTomato] with the following enzymes: NheI-HF and Ascl.

Plasmid 84484:
pAAV-CAG-FLEX-rc [Chronos-tdTomato]



DNA digested with NheI-HF and Ascl single digests and NheI-HF/ Ascl double digest for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.