


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pEN395 - pCAGGS-Tir1-V5-BpA-Frt-PGK-EM7-PuroR-bpA-Frt TIGRE targeting with the following enzymes: EcoRV-HF/ HindIII-HF

Plasmid 92141

pEN395 - pCAGGS-Tir1-V5-BpA-Frt-PGK-EM7-PuroR-bpA-Frt
TIGRE targeting



DNA digested with EcoRV-HF/ HindIII-HF double digests for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.